

XAFS Tutorials for crystallographers and beginners: General introduction and analytical methods for crystallographers and beginners in the field of XAFS.

FEFF and Related Codes

Anatoly Frenkel

Professor

Physics Department, Yeshiva University, New York, USA

Synchrotron Catalysis Consortium, Brookhaven National Laboratory, USA

www.yu.edu/faculty/afrenkel

Anatoly.Frenkel@yu.edu

FEFF: John Rehr (U. Washington) and his group members


IFEFFIT: Matthew Newville (U. Chicago) and Bruce Ravel (NIST)

AUTOBK: Matthew Newville


(extension of the original UWXAFS package developed
under the direction of Edward Stern)

Capabilities of EXAFS for studies of nanoparticle catalysts


-Size:



-Structure:



-Shape:




- Short range order:

-Texture:



-Core/shell or surface segregation




A. I. Frenkel, A. Yevick, C. Cooper, R. Vasic
Ann. Rev. Anal. Chem., **4**, 23-39 (2011)

A. I. Frenkel, Z. Kristallographie, **222**, 605-611 (2007).

A. I. Frenkel, C. W. Hills, and R. G. Nuzzo,
Feature Article, J. Phys. Chem. B, **105**, 12689-12703
(2001).

A. I. Frenkel, J. Synchrotron Rad., **6**, 293 (1999).

-Random nanoalloys



Theoretical EXAFS Equation

Single scattering path:

$$\chi_{\Gamma}(k) = \frac{NS_0^2}{kR^2} \left| f^{\text{eff}}(k) \right| e^{-2\sigma^2 k^2} e^{\frac{-2R}{\lambda}} \sin \left[2kR - \frac{4}{3} C_3 k^3 + \delta(k) \right]$$

Multiple-scattering path:

$$\chi_{\Gamma}(k) = \text{Im} \, NS_0^2 \frac{e^{i \left(k \sum_i R_{ii+} + 2\delta(k) \right)}}{\prod_i kR_{ii+}} e^{-2\sigma^2 k^2} e^{\frac{-2R}{\lambda}} \text{Tr} \, MF^N \dots F^2 F^1$$

Path expansion: $\chi(k) = \sum_{\Gamma} \chi_{\Gamma}(k)$

Parameters in EXAFS Equation

$$\chi_{\Gamma}(k) = \frac{N S_0^2}{k R^2} |f^{\text{eff}}(k)| e^{-2\sigma^2 k^2} e^{\frac{-2R}{\lambda}} \sin\left[2kR - \frac{4}{3} C_3 k^3 + \delta(k)\right]$$

Amplitude → $N S_0^2$
FEFF → $f^{\text{eff}}(k)$
Amplitude → σ^2
FEFF → λ
Amplitude → C_3
FEFF → $\delta(k)$

Same parameterization for SS and MS paths.

$$\chi(k) = \sum_{\Gamma} \chi_{\Gamma}(k)$$

$$k = \sqrt{\frac{2m}{\hbar^2} (E - E_0)}$$

$$E_0 = E_0^{\text{bkg}} + \Delta E_0$$

$$R = R_{\text{model}} + \Delta R$$

IFEFFIT
 (GUI: Artemis)

- Easy and transparent parameterization of structural models for EXAFS data fitting
- FEFF theory
- Error analysis

Main similarities and differences between EXAFS analysis packages

Main purpose:

- Construct structural model,
- Calculate theoretical EXAFS signal,
- Perform structural refinement,
- Carry out adequate error analysis,
- Pick the model that fits the data the best and **use the best fit results**
- Keep in mind that EXAFS is ensemble averaging method!

Use complementary techniques

Difference between packages:

Fitting Space: k-space or r-space

Parameterization of a model

Multiple data set fitting

Theoretical method
(muffin tin or finite difference method)

FEFF

FEFF6,8, 9:

Spherically symmetric potentials (muffin tin approximation)

Input: (xyz) coordinates and atomic numbers

Mean free path

Imaginary part of interstitial potential and life time broadening

Self Energy

Metals – Hedin-Lundqvist

Insulators – Hedin-Lundqvist or Dirac-Hara

Molecules – Dirac-Hara or ground state

Multiple scattering expansion

Each photoelectron path with its $f(k)$, $\delta(k)$, $\lambda(k)$ is saved as a file

Packages using FEFF:

IFEFFIT, EXAFSPAK, WINXAS, XDAP,

XFIT, EDA, LASE, MAC

Fitting EXAFS Theory to the Data:

$$f(R_i) = \tilde{\chi}(R_i) - \tilde{\chi}_M(R_i)$$


$$\chi_\nu^2 = \frac{1}{\nu} \sum_{i=1}^{N_{\text{idp}}} \left(\frac{f_i}{\varepsilon_i} \right)^2 = \frac{N_{\text{idp}}}{N\nu\varepsilon^2} \sum_{i=1}^N \left[\text{Re}(f_i)^2 + \text{Im}(f_i)^2 \right]$$

$$\nu = N_{\text{idp}} - P \quad (\text{Number of degrees of freedom})$$

$$N_{\text{idp}} = \frac{2\Delta k \Delta R}{\pi} \quad (\text{Number of relevant independent data points})$$



EXAFS data analysis and modeling

I. Processing and visual examination of the data




II. Deciding on the model and refinement parameters


Structural modeling
(different for homogeneous and heterogeneous systems)



Path expansion and FEFF modeling
 $f(k)$, $\delta(k)$, $\lambda(k)$



III. Fitting theoretical EXAFS signal to the data



Step by step data analysis tutorial

I. Analysis of bulk Pd foil

Directory: Foil fitting

Copy the Pdfoil.apj into Pdfoil-SingleShell.apj

Start Artemis application

Open Pdfoil-SingleShell.apj into Artemis

Discard existing theoretical signal (we will construct a new one from scratch) by

Right click on FEFF0, follow the link to theory → discard this FEFF calculation.

Click on Theory menu → New Atoms page → Blank page

Space group: f m 3 m (remember to put space between symbols)

A: 3.89

Element: Pd

Tag: Pd

X: 0

Y: 0

Z: 0

Click Define

Click Run Atoms

Click Run Feff

Pick “Just the first” option (the path1 will appear in the central panel).

Highlight the Path1 and you will see options in the left panel. Give names S02, enot, drpd, ss to the variables amp, enot, delr, ss, respectively (or any other names)

Press Guess Def Set

Delete all values that are pre-set there

In the central panel, highlight the Path1 again

Right click on each variable and pick “make a guess and stay”

Go back to Guess Def Set to make sure all the variables appear there. In the S02 value change the guess value from 0 to 1 and press enter

1. Analysis of bulk Pd foil (Continued)

Highlight the data and press R in the right panel

Adjust the r-range for the fit (make r-range from 1 to 3.06) to define the region of the first peak

Press the Fit button

Examine the fit quality in the graphic window. If it is good examine the results in the Palette.

Take a note of the best fit value for S02: 0.825

2. Analysis of Pd nanoparticle data

Open project file Pd-nw.....

Click on the Guess Def Set and replace the pre-set value of 0.815 by the value we obtained in the fit to the foil: 0.825

The theory now has the amplitude factor written in the form: the product of S02 (which we found from the fit to the foil) and the unknown coordination number Npd (which was fixed at 12 when we analyzed the foil). The rest of the fitting parameters are the same.

Press Fit and examine the fit quality and the results.

The coordination number is 9.31 ± 0.35 . That corresponds to a particle size of around 1.7 nm in diameter (assuming cuboctahedral geometry and particle on support).

Questions: contact anatoly.frenkel@yu.edu