

IUCr-XAFS, Madrid, 22-8-2011

APPLICATIONS OF X-RAY SPECTROSCOPY IN CHEMISTRY

M. Fernández-García

Instituto de Catálisis y Petroleoquímica
ICP-CSIC, Spain

APPLICATIONS OF X-RAY SPECTROSCOPY IN CHEMISTRY

- **Summary**

- General Introduction XAFS/XES

Overview XANES/EXAFS utility in Chemistry

- XAFS analysis of materials genesis

- XAFS analysis of industrial-related chemical properties

APPLICATIONS OF X-RAY SPECTROSCOPY IN CHEMISTRY

- **XAFS (XES) Techniques: XANES and EXAFS**

Fit Chemical needs (Atmosphere, P, T)

Local order sensitive (around each component)

Independent of crystallinity

XANES: high S/N ratio (absence thermal effects)

NANO-CHEMISTRY

Control of Structure and Morphology (primary particle size, shape)


Structure / Morphology \leftrightarrow Chemical Properties

OVERVIEW XANES APPLICATION: CHEMISTRY

Electronic Properties

- Edge
- Pre-edge
- CRs

$$\Delta E = \Delta \sum_j^{\text{val}} \left\langle e^2 / r_j \right\rangle + E(\text{exchange})$$




Average Oxidation State

OVERVIEW XANES APPLICATION: CHEMISTRY

Electronic Properties

- Edge
- Pre-edge
- CRs

W: TiO₂ W L_I L_{III}-edge



Pre-edge Position and Intensity: Ox. State + Local symmetry


CRs: “Projected” Empty (s,p,d) Density of States

OVERVIEW XANES APPLICATION: CHEMISTRY

Electronic/Structural Properties

- Pre-edge
- CRs
- Nanophases
- Th. Calculations
- Size/Shape dependence


Cu nanoparticles




OVERVIEW XANES APPLICATION: CHEMISTRY

Electronic Properties

- Edge
- Pre-edge
- CRs



W: TiO_2 XES W L_{III}-edge




High Resolution => Fine electronic details

OVERVIEW XAFS APPLICATION: CHEMISTRY

XES: selective spectra

Electronic Properties


- XANES




OVERVIEW XAFS APPLICATION: CHEMISTRY

Structural Properties

- EXAFS



XES: selective spectra



Site (Local St.) Selective Structural details

OVERVIEW XANES APPLICATION: CHEMISTRY

Structural Properties


- Pre-edge ($s \rightarrow p/d$)

Undercoordinated
Metal



Number of Oxygen
Vacancies

Metal Oxides (1st Transition Series)



OVERVIEW XANES APPLICATION: CHEMISTRY

Metal Oxides (1st Transition Series)

Structural Properties


- Pre-edge (s \rightarrow p/d)

Partially Reduced
And/or
Undercoordinated
Metal



Number
Defects/Vacancies

X-Ray Spectrom. 37 (2008) 572




OVERVIEW XANES APPLICATION: CHEMISTRY


Structural Properties

Ce^{3+} in quasi-stoichiometric particles

- CRs



- $\text{Ce}^{3+} \leftrightarrow$ Vacancies below 5 at. %



- Strain from Ce^{3+} presence
Raman ($\text{Ce}^{3+} < 1.6$ at. %)

XANES: Relatively Poor Detection Limit

PCCP 6 (2004) 3524

OVERVIEW XANES APPLICATION: CHEMISTRY

Structural Properties

- CRs

$$E \propto 1/R^2$$

Phase behavior



Analysis
CR 9020 eV



disordered FCC
alloys

Local ordering
 $Pd/Cu \leq 1$

Pd-Cu alloys



OVERVIEW EXAFS APPLICATION: CHEMISTRY

Structural Properties

- Primary p. size and shape for nano-particles



1st Shell CN

Size

Dispersion (Surf. vs. Bulk)

High Sensitivity < 300 at.

Nano-metals



OVERVIEW EXAFS APPLICATION: CHEMISTRY

Structural Properties

- Primary p. size/shape

1st Shell CN

Size


Surface Disorder

Size < 5 nm

Limited validity
Usual Analytical Approach

Underestimation CN

Nano-metals



OVERVIEW EXAFS APPLICATION: CHEMISTRY

Structural Properties

- Primary p. size and shape for nano-particles (metals)


1st Shell CN; size

CN_n: CN₁; shape

Fcc; hcp OK

Bcc; only CN₁ > 4

Nano-metals



OVERVIEW EXAFS APPLICATION: CHEMISTRY

Structural Properties

Homogeneity

$$\alpha = \frac{1 - N_{\text{PtPd}} / N_{\sum \text{PtPt} + \text{PtPd}}}{X_{\text{Pd}}}$$

$$\Delta = \alpha_{\text{M1}} - \alpha_{\text{M2}} \quad \Sigma = \alpha_{\text{M1}} + \alpha_{\text{M2}}$$

Combination

Degree of interaction/mixing

Nano-alloys




Fig. 9 Hypothetical bimetallic clusters (Pt = blue, Pd = green): (a) bi-cluster; (b) core-shell; (c) random; (d) single-half; (e) bi-cluster homo; (f) bi-cluster inhomogeneous.

Table 5 1st shell coordination numbers for PtPd bimetallic systems

Cluster Type	Pt-Pt	Pt-Pd	$\sum N_{\text{PtP}}$	Pd-Pd	Pd-Pt	$\sum N_{\text{PdP}}$	x_{Pt}	x_{Pd}	α_{Pt}	α_{Pd}	$\Delta\alpha$	$\sum\alpha$
Bi-cluster	8.73	0.07	8.80	8.73	0.07	8.80	0.5	0.5	-0.98	-0.98	0	-1.96
Single-Half	7.08	1.58	8.66	7.85	1.19	9.04	0.43	0.57	-0.68	-0.69	0.01	-1.37
Core-shell	7.73	3.86	11.59	3.90	2.92	6.82	0.43	0.57	-0.42	0	-0.42	-0.42
Random	3.63	4.80	8.42	5.59	3.63	9.22	0.43	0.57	0	-0.08	0.08	-0.08
Decorated	3.43	3.00	6.43	8.73	0.77	9.50	0.21	0.78	-0.18	-0.9	0.72	-1.08
Random (50/50)	5.41	4.26	9.68	4.20	3.88	8.09	0.5	0.5	-0.12	-0.08	-0.04	-0.2
Bi-cluster homo	7.63	0.92	8.55	3.25	4.33	7.58	0.82	0.18	-0.39	-0.31	-0.08	-0.7
Bi-cluster inhomogeneous	7.63	0.24	7.88	5.92	1.17	7.08	0.82	0.18	-0.82	-0.8	-0.02	-1.62

APPLICATIONS OF X-RAY SPECTROSCOPY IN CHEMISTRY

- Summary

- General Introduction XAS/XES

Overview XANES/EXAFS utility in Chemistry



- *XAFS analysis of materials genesis*

- **XAFS analysis of industrial-related chemical properties**

CONTROL OF STRUCTURE AND MORPHOLOGY: SYNTHESIS

Liquid Phase; synthesis of Zeotypes (Co-AlPO)

XAS - XRD Analysis




Complementary analysis of the local/long range order evolution

CONTROL OF STRUCTURE AND MORPHOLOGY: SYNTHESIS

Liquid Phase; synthesis of Zeotypes

XAS XRD Analysis



Detection of four different regions
(Initial) Intermediate only detected by XAS


CONTROL OF STRUCTURE AND MORPHOLOGY: SYNTHESIS

Synthesis of Oxides

Surfactant	Triton X-100			Tween 85	
H ₂ O / Surf.	18	4.5	3	18	3
Samples / 773 K	T	TA	TB	TW	TWB
Samples / 873 K	T1	T1A	T1B	T1W	T1WB

Solid state treatment


High Surface Area
Nanostructured Anatase materials



CONTROL OF STRUCTURE AND MORPHOLOGY: SYNTHESIS

Solid phase: Initial Precursors


XANES Spectra: Pre-edge



Anatase-type Structure



EXAFS Spectra




Ti-Ti C.N. 2.2 ± 0.5

Precursors lack 3D connectivity

CONTROL OF STRUCTURE AND MORPHOLOGY: SYNTHESIS

Treatment of Solid precursors; synthesis of Oxides

EXAFS Spectra: Ti-Ti 2nd Shell CN





$R < 4 \text{ \AA}$ may inform about the *nucleation* step

CONTROL OF STRUCTURE AND MORPHOLOGY: SYNTHESIS

Treatment of Solid precursors; synthesis of Oxides

Multitechnique XAS; XRD-PDF



3-5 Å local order → control *nucleation* and primary particle size

Thermodynamic factor


CONTROL OF STRUCTURE AND MORPHOLOGY: SYNTHESIS

Treatment of Solid precursors; synthesis of Oxides

XAS Analysis

$R > 4 \text{ \AA}$ may inform about the *growth* process (at low T)


XRD Analysis



Growth ($D^2 - D_0^2$) = kT ; 15-30 kJ mol^{-1} ; hydration surface layers

CONTROL OF STRUCTURE AND MORPHOLOGY: SYNTHESIS


Surface “sensitive” DRIFT Spectra




Control particle shape: complex phenomenon (thermo/structure)

CONTROL OF STRUCTURE AND MORPHOLOGY: SYNTHESIS

TiO₂ Anatase Nanomaterials: Size and Shape





Shape



Morphology control :

Size (Triton) Shape (“Tween”)

CONTROL OF STRUCTURE AND MORPHOLOGY: ACTIVITY



3 Size Regions

Size < 5-6 nm; low activity

NO clear differences a/b 15 nm

Shape

Above 6 nm

Isotropic Shape enhances Photoactivity

APPLICATIONS OF X-RAY SPECTROSCOPY IN CHEMISTRY

- Summary

- General Introduction XAFS/XES


Overview XANES/EXAFS utility in Chemistry

- XAFS analysis of materials genesis

- *XAFS analysis of industrial-related chemical properties*

CHEMICAL ACTIVITY OF MATERIALS

Mutitechnique SYNCHROTON approach



XAFS-SAXS-WAXS; XAFS-UV (MS) approaches

Real “operando” conditions





Time dependent Structure – Chemical Activity

Catal. Today 145 (2009) 204

CHEMICAL ACTIVITY OF MATERIALS

Mutitechnique SYNCHRONOUS approach Modern requirements




Energy Dispersive ED-XAFS-IR (MS)

μs “real” single-shot Time Resolved experiment; TOFs

Catal. Today 145 (2009) 204

CHEMICAL ACTIVITY OF MATERIALS

Mutitechnique SYNCHRONOUS approach Modern requirements




Spatial Resolved ED-XAFS-IR (MS)

2D (3D) sub- μm analysis of industrial Structures
(pellets, monolith)

Chem. Soc. Rev. 39 (2010) 4741

CHEMICAL ACTIVITY OF MATERIALS XAS techniques; novel experiment



Summary



Fast (ms- μ s) easy track of structural-electronic properties

Combination Surface-Bulk
(Spatial-resolved) Information

Link with Chemical Activity

CHEMICAL ACTIVITY OF MATERIALS

XANES technique; time-resolved studies




Size-dependence => Fitting procedures

Statistical Analysis: Self-resolving or multivariate curve resolution

Angew. Chem. 46 (2007) 8629

CHEMICAL ACTIVITY OF MATERIALS

XANES-Fluorescence microscopy technique; Spatial Resolution




2D Chemical (10 x 10 nm x 10-100 μm) imaging of
“active” species

Nature 456 (2008) 222

CHEMICAL ACTIVITY OF MATERIALS

XANES (XRD) microtomography technique; Spatial Resolution



3D Chemical (10x10 nm x 360 °) imaging of “active” species


Appl. Phys. Lett. 82 (2003) 3360

CHEMICAL ACTIVITY OF MATERIALS EXAFS technique; time-resolved studies

Pd size-dependence; cycling CO/NO at 673 K

EXAFS analysis assumptions



σ^2 constant and/or R harmonic distribution




Incorrect Information

J. Am. Chem. Soc. 132 (2010) 4540

CHEMICAL ACTIVITY OF MATERIALS
EXAFS technique; time-resolved studies
Pd size-dependence; cycling CO/NO at 673 K




XAS; XRD Information



J. Am. Chem. Soc. 132 (2010) 4540

J. Catal. 270 (2010) 275

CHEMICAL ACTIVITY OF MATERIALS
EXAFS technique; time-resolved studies
Pd size-dependence; cycling CO/NO at 673 K




CHEMICAL ACTIVITY OF MATERIALS


EXAFS technique; time-resolved studies

Pd size-dependence; cycling CO/NO at 673 K

Raw Data



Dynamic Dispersion



“Dynamic” Surface Normalized

(TOFs) “True”Activity

APPLICATIONS OF X-RAY SPECTROSCOPY IN CHEMISTRY

XANES/EXAFS

Tools for dynamic structural-electronic characterization of
nano-solids

CONTROL PHASE AND MORPHOLOGY

CHEMICAL ACTIVITY

ACKNOWLEDGEMENTS

Collaborators

Dr. A. Kubacka
Dr. Iglesias Juez
Dr. K. Chistoforidis
Dr. M. Ferrer
Prof. A. Martínez Arias
Prof. J.C. Conesa
Prof. J. Soria
ICP; CSIC; SPAIN

Dr. G. Colón
Prof. J.A. Navío
ICMSE; CSIC-US; SPAIN

Dr. J.A. Rodríguez
Dr. J.C. Hanson
BNL; USA

Dr. S. Figueroa
Dr. M.A. Newton
Dr. M. Di Michiel
ESRF; France

XAS experiments

All staff at:
-7.1, 9.2, 9.3 at SRS

*-ID15, BM23, ID24, ID26,
BM29 at ESRF*

-X7B, X19A at NSLS

Financial sources

EU program for Large Scale Installations
Spain-MICINN (CTQ2010-14872)